



## DDS Full Comparison Review

Brand	Easy-to-Read Gauge	Patents	Cleanliness	Sizes	Lumbar Curvature Design
<b>DDS</b>	<b>YES</b>	<b>5</b>	Machine washable Hand wash Air dry in the shade	<b>6 +</b>	<b>YES</b>
<b>Dr. Ho</b>	<b>NO</b>	<b>1</b>	Clean with warm water and soap. Air dry only. Lay flat to dry. Do not submerge in water. Do <b>NOT</b> use dryer. Do <b>NOT</b> iron. Do <b>NOT</b> bleach. Do <b>NOT</b> put in direct sunlight.	<b>2</b>	<b>NO</b>
<b>Alphay</b>	<b>YES</b>	<b>1</b>	<b>NO</b>	<b>2</b>	<b>NO</b>

Brand	Easy-to-Read Gauge	Patents	Cleanliness	Sizes	Lumbar Curvature Design
HailiCare	NO	1	<p>Don't immerse into water to wash it.</p> <p>Avoid direct sunlight and heat.</p> <p>Put it in a cool, dry place after using.</p>	1	NO
EHM	YES	0	No water	0	NO
No Brand	NO	0	No water	2	NO
Tcare Physio	YES	0	No water	2	NO

**Easy-to-Read Gauge:** Provides an accurate reading for personalized and consistent treatment each and every time. Individuals have the option to use customized levels of pressure during specific activities such as driving or during heavy lifting; or during different times of the day.

**Lumbar Curvature Design:** The design allows the belt to conform to the natural curvature of the lumbar spine as it expands vertically whereas belts without this feature inflates flat and straight. The designed curve is essential to effectively anchor underneath the lowest rib pushing upwards and against the hip pushing down during expansion. This action gently stretches the torso vertically and displaces pressure away from the affected disc and nerve. Flat or straight expansion may not properly anchor on these areas, therefore impeding on its' efficacy.